

# SONIC


DETACHABLE TIP BRAIDED MICROCATHETER  
COMPATIBLE WITH DMSO



## DETACHABLE TIP

The patented FUSECATH system will allow catheter tip to be detached when trapped, enabling the catheter to be removed.

## FUSECATH® SYSTEM




### 3 radiopaque markers

- 1 • Distal tip
- 2 • Detachment point
- 3 • Maximum reflux



# SONIC

## FLOW-DIRECTED BRAIDED MICROCATHETER




## SONIC RANGE

REFERENCE	TOTAL LENGTH	DISTAL TIP O.D.	DISTAL TIP I.D.	FUSECATH LENGTH	DEAD SPACE	MAX REFLUX	INCLUDED GUIDEWIRE
SONIC1.2F15	165cm	1.2F	0,17mm	1,5cm	0,3ml	2cm	.007"
SONIC1.2F25	165cm	1.2F	0,17mm	2,5cm	0,3ml	3cm	.007"
<b>NEW!</b> SONIC1.2F35	165cm	1.2F	0,17mm	<b>3,5cm</b>	0,3ml	4cm	.007"
SONIC1.2F25.190	190cm	1.2F	0,17mm	2,5cm	0,35ml	3cm	.007"
<b>NEW!</b> SONIC1.5F15	165cm	1.5F	0,25mm	<b>1,5cm</b>	0,33ml	2cm	.008"
SONIC1.5F25	165cm	1.5F	0,25mm	2,5cm	0,33ml	3cm	.008"
<b>NEW!</b> SONIC1.5F35	165cm	1.5F	0,25mm	<b>3,5cm</b>	0,33ml	4cm	.008"

All SONIC are provided with a HYBRID guidewire.

## TRI-AXIAL APPROACH

For a tri-axial approach, use SONIC catheter together with FARGOMAX and FARGOMINI.



### CAUTION

SONIC catheter should be used only with the indicated guidewire.



10, rue de la Croix Vigneron  
95160 Montmorency • France  
Tél. : +33 (0)1 39 89 46 41  
Fax : +33 (0)1 34 17 03 46  
www.balt.fr